

Thopaz™



## PATIENT TUBING OPTIONS




Thoracic drainage system







# Thopaz™ Tubing

Thopaz double lumen tubing is available in single or double patient connector configurations that accommodate a variety of chest catheter sizes. The tubing helps measure and regulate prescribed negative pressure close to a patient's chest while removing chest secretions for collection in an attached canister.

In addition, Thopaz monitors the tubing to help prevent clogging and contains a sample port for convenient specimen collection of drainage fluid. (Note: For sample collection, use a 17G or thinner blunt cannula).



## Thopaz Tubing

	0790023	Small Thopaz Tubing, Single Connector with Sample Port. Compatible with 12-20Fr catheter.	10/CS
	0790024	Small Thopaz Tubing, Double, Connector with Sample Port. Compatible with 12-20Fr catheter.	10/CS
	0790021	Medium Thopaz Tubing, Single Connector with Sample Port. Compatible with 20-32Fr catheter.	10/CS
	0790022	Medium Thopaz Tubing, Double Connector with Sample Port. Compatible with 20-32Fr catheter.	10/CS
	0790025	Large Thopaz Tubing, Single Connector with Sample Port. Compatible with 32-40Fr catheter.	10/CS
	0790026	Large Thopaz Tubing, Double Connector with Sample Port. Compatible with 32-40Fr catheter.	10/CS

Medela offers healthcare professionals 50 years of expertise in medical vacuum technology. Our Swiss development and manufacturing process yields mechanical precision to provide excellence in performance.

Medical Vacuum Technology for Healthcare Professionals

Please contact us or your local Medela representative for details.

**CAUTION: Federal US law restricts this device to sale by or on the order of a physician.**

 Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
www.medela.com

Medela Inc. USA  
1101 Corporate Drive  
McHenry, IL 60050  
USA  
Phone +1 877 735 1626  
Fax +1 815 307 8942  
suction@medela.com  
www.medelasuction.com

Medela is a registered trademarks of Medela Holding AG.  
Thopaz is a trademark of Medela, Inc.

Local contact: